ZIMMER BIOMET

X Series Power System Owner's Manual

TABLE OF CONTENTS

GENERAL INFORMATION	4
OWNER'S MANUAL SCOPE	4
INTENDED USE	4
CONTRAINDICATIONS	4
GENERAL CAUTIONS	4
INITIAL VISUAL INSPECTION	5
HANDPIECES	6
FUNCTIONALITY CHECK	6
ATTACHMENTS	9
ATTACHMENT INSERTION AND DISCONNECTION	9
ROTARY ATTACHMENTS	9
HIGH TORQUE REAMER ATTACHMENTS	9
STANDARD TORQUE REAMER ATTACHMENTS	10
DRILLING ATTACHMENTS	10
SAW ATTACHMENTS	13
SAW ATTACHMENTS	13
SPECIALTY ATTACHMENTS	16
STERNUM SAW GUARDS	16
PIN AND WIRE DRIVERS	17
POWER SOURCES AND CHARGERS	19
BATTERIES	19
STERILIZABLE BATTERY	19
ASEPTIC BATTERY	19
SETUP	20
BATTERY CHARGERS	20
CHARGER FOR ASEPTIC BATTERY	20
CHARGER FOR STERILIZABLE BATTERY	20
SETUP	21
FUNCTIONALITY CHECK	23
STERILIZATION BASKETS	25
STERILIZATION BASKET FOR USE WITH ASEPTIC BATTERY HOUSINGS	25
SYSTEM FUNCTIONS	26
OPERATION	26
SURGEON DECISIONS	26
TESTS AND CHECKS	26
MAINTENANCE INDICATORS	26
ASSEMBLY OF THE HANDPIECE, ATTACHMENT AND BATTERIES	27
ASEPTIC BATTERY TRANSFER	28
CLEANING, DECONTAMINATION AND STERILIZATION	29
PRECAUTIONS FOR CLEANING, DECONTAMINATION AND STERILIZATION	30
PREPARING FOR MANUAL CLEANING	30
MANUAL CLEANING MATERIALS	30
MANUAL CLEANING PROCESS	31
MANUAL CLEANING INSPECTION	31
DISINFECTION	31
PREPARING FOR AUTOMATED CLEANING	31
AUTOMATED CLEANING MATERIALS	32
AUTOMATED CLEANING PROCESS	32
AUTOMATED CLEANING INSPECTION	32

PREPARING FOR STERILIZATION	32
STERILIZATION OF HANDPIECE, ATTACHMENTS, ASEPTIC BATTERY HOUSING AND ASEPTIC BATTERY FUNNEL	33
STEAM STERILIZATION	33
IMMEDIATE USE STEAM STERILIZATION (IUSS)	33
STERILIZATION OF STERILIZABLE BATTERY	34
STEAM STERILIZATION	34
IMMEDIATE USE STEAM STERILIZATION (IUSS)	34
STERRAD	34
V-PRO	34
MAINTENANCE	35
GENERAL MAINTENANCE INFORMATION	35
SERVICING	35
PRECAUTIONS IN SERVICE	35
PROTECTIVE EARTH RESISTANCE TESTING	35
EARTH (EQUIPMENT) LEAKAGE CURRENT TESTING	35
INSULATION RESISTANCE TESTING	36
BATTERY SERVICE	36
BATTERY-HANDLING RECOMMENDATIONS	36
FUSE AND FUSE DRAWER REPLACEMENT	37
UNSCHEDULED MAINTENANCE	37
TROUBLESHOOTING GUIDE	37
TROUBLESHOOTING	37
EXPECTED SERVICE LIFE INFORMATION	38
RETURN AUTHORIZATION AND REPLACEMENT INFORMATION	38
REPLACEMENT PARTS	38
WARRANTIES	39
WARRANTY INSIDE THE U.S.A.	39
DISCLAIMER OF IMPLIED WARRANTIES	39
LIMITATION OF REMEDIES	39
WARRANTY CLAIMS	39
WARRANTY OUTSIDE U.S.A.	39
SCOPE OF WARRANTY	39
SPECIFICATIONS AND PERFORMANCE	40
CHARACTERISTICS	40
IEC 60601-1 Classification	41
EMISSIONS/IMMUNITY	41
ELECTROMAGNETIC COMPATIBILITY (EMC) GUIDANCE TABLES	41
EMC GUIDANCE AND DECLARATION – EM EMISSIONS	42
EMC GUIDANCE AND DECLARATION – EM IMMUNITY/DISTURBANCES	42
EMC GUIDANCE AND DECLARATION – EM EMISSIONS/RF	43
EMC GUIDANCE AND DECLARATION – IMMUNITY/SEPARATION DISTANCES	43
ENVIRONMENTAL CONDITIONS	44
SYMBOLGY	45
TEST AND EVALUATION REPORT	47
CONTACT INFORMATION	48
REPORTING PROBLEMS	48

GENERAL INFORMATION

OWNER'S MANUAL SCOPE

This Owner's Manual specifies characteristics related to the use of the X Series Power System as well as its intended use, intended patient population as prescribed by the clinician, intended user profile (certified clinician), intended conditions of use as defined in this manual and operating principles. The scope of this manual will address the handpieces, the attachments, the batteries, the battery chargers and the sterilization baskets.

In addition to this Owner's Manual, each product is also provided with a product-specific instruction manual contained within each package.

INTENDED USE

The intended use of the X Series Power System is for cutting, shaping, and preparing bone and bone-related tissues in a variety of surgical procedures. The system is usable for performing numerous operations such as sawing, reaming and drilling as well as driving wires, pins and various instruments.

CONTRAINDICATIONS

None known.

CAUTION: In every case, the final decision whether to use the X Series Power System rests with the attending physician. He/she should assess prior to use of the system if its performance specifications, precision and handling are suitable for a given application, taking into account the bone strength, anatomical situation and bone size. Moreover, the technique guides and any contraindications of the implant system must be respected.

GENERAL CAUTIONS

- Only medical professionals who are thoroughly familiar with this power equipment's function, application and instruction for use should operate.
- Do not use if the device appears damaged or does not work as intended.
- Cutting tools must be discarded in sharps container after use.
- Not for use in an oxygen rich environment with an oxygen concentration greater than 25% for ambient pressures up to 110 kPa or where the partial pressure of oxygen is greater than 27.5 kPa at ambient pressures exceeding 110 kPa.
- Do not use accessories or cords other than those specified or sold by Zimmer Biomet.
- Use only with attachments and power sources that are part of the X Series Power System.
- Use hospital grade plug with the minimum specifications listed in the specifications section of this manual.
- Grounding reliability can only be achieved when the equipment is connected to an equivalent receptacle marked hospital only or hospital grade.
- Do not use the battery charger if it has become saturated or immersed in fluid.
- Ensure the power cord does not create a trip hazard.
- Never block the air vents of the battery charger.

GENERAL INFORMATION

- Do not immerse the battery charger or batteries in fluids.
- The ON/OFF switch on the battery charger will not completely shut off all power from the AC outlet. Since the power cord serves as the main disconnect device for the unit, unplug the power cord from the AC outlet to shut down all power. Make sure the unit has been installed so that the main power cord can be easily unplugged from the AC outlet if needed.
- It is highly recommended to continuously moisten bones when drilling, cutting or reaming.
- When the handpiece is not in use, place the mode switch in the center position to lock-out the trigger(s) and remove the power source.
- In the event the handpiece is left connected to a power source for a period longer than 30 minutes, the system will enter a low power mode. To return to normal mode, pull the lower trigger once (for a double trigger handpiece) or the single trigger once (for a single trigger handpiece) and the system will return to normal mode. In the event the handpiece is left connected to a battery for a period longer than 45 minutes, the battery will enter low power mode and the battery LEDs will turn OFF. To return to normal mode, pull the lower trigger once (for a double trigger handpiece) or the single trigger once (for a single trigger handpiece) and the battery will return to normal mode.
- Coupling and decoupling the attachments from the handpiece or coupling and decoupling the adaptors or cutting tools (drill bits, reamers, saw blades) from the attachments must not be performed above the patient.
- Always manipulate or change the attachments and accessories on the back table with the handpiece mode switch set to the center locked position.
- Use only saw blades sold by Zimmer Biomet.
- To the extent possible, avoid contact with the cutting tools (drill bits, reamers, saw blades) and adapters once the handpiece is connected to a power source.
- Do not touch the terminals of the batteries with metallic objects.
- Do not use alkaline cleaning agents with a pH value greater than 10.8.
- For cleaning and decontamination, use only detergents compatible with anodized aluminum surfaces.
- The power equipment should never be immersed in liquid.
- The power equipment should never be cleaned in an ultrasonic cleaning system.
- Do not place the batteries or battery chargers in an automated washer.
- Do not sterilize the red aseptic battery or battery chargers.
- Exposing the aseptic battery to extreme temperatures may reduce battery capacity and/or performance. Extreme temperatures could result in a fire or burn hazard.
- The empty aseptic battery housing must be washed and sterilized with its lid in an open position.
- Allow the handpiece, the sterilizable battery, the aseptic battery housing and attachments to properly cool after removing from the autoclave to normal room temperature before using or storing.
- A heated handpiece can be used normally but should not be left in continued contact with the patient's tissues when not used for its intended purpose.
- The handpieces, batteries and battery chargers contain electrical materials. The presence of these materials may, if not disposed of properly, have potential adverse effects on the environment. Do not dispose of in normal household waste.
- The manual cleaning and disinfection procedure must be performed prior to the automated cleaning and disinfection procedure. Use a non-abrasive, non-metallic brush to clean all surfaces including the complex surfaces before the automated wash cycle.
- Prior to initial usage, the X Series Power System must be cleaned and sterilized per the instructions (See Cleaning, Decontamination and Sterilization Instructions).
- The use of Zimmer Biomet cleaning and sterilization baskets is highly recommended.

INITIAL VISUAL INSPECTION

Unpack the X Series Power System upon receipt and inspect the unit for any obvious damage that may have occurred during shipment including damage to any part of the system. It is recommended that this inspection be performed by a qualified biomedical engineer or other person thoroughly familiar with electronic medical devices. If the unit is damaged, notify the carrier and your Zimmer Biomet representative immediately.

HANDPIECES

REF 89-8520-400-00

REF 89-8520-400-10

The X Series Power System handpieces have a variety of features and physical characteristics as described below.

Double Trigger Handpiece

Single Trigger Handpiece

	Feature	Action
1	Main Body	Main body of the handpiece that holds the electronics and precision drive motor.
2	Attachment Release	Attachment release mechanism that, when pushed, allows the attachment to be removed.
3	Mode Switch	Mode switch that allows both the double trigger and single trigger systems to have the triggers locked in the OFF position. For the DOUBLE TRIGGER handpiece, the mode switch operates as: FORWARD-OFF-FORWARD For the SINGLE TRIGGER handpiece, the mode switch operates as: REVERSE-OFF-FORWARD
4	Double Trigger UPPER TRIGGER	Trigger for providing the REVERSE (counterclockwise) function of the double trigger handpiece. When this trigger is held down along with the LOWER trigger, this activates the OSCILLATION feature.
5	Double Trigger LOWER TRIGGER	Trigger for providing the FORWARD (clockwise) function of the double trigger handpiece. When this trigger is held down along with the UPPER trigger, this activates the OSCILLATION feature.
6	Single Trigger TRIGGER	Trigger for providing the FORWARD and REVERSE function of the single trigger handpiece. The mode switch is used to change the rotational direction.

CAUTION: Ensure familiarity with the above trigger functionality prior to use.

FUNCTIONALITY CHECK

This section describes automatic checks and manual checks to verify the handpiece functions. The automatic diagnostics are built into the X Series Power System and will check the select software and hardware functions upon power-up. This automatic check will take place when the power is connected. A passing result for the automatic checks allows the handpiece to drive the motor forward, reverse and oscillate modes; a failing result stops the handpiece from working. If the automatic checks pass, the manual checks can be performed. The manual checks can be performed by following the instructions below.

CAUTION: It is recommended to verify the functionality of the handpiece prior to each use.

WARNING: To avoid the risk of electric shock, the battery chargers must only be connected to an ~ AC Power Mains with protective earth.

1. Verify the handpiece MODE SWITCH is set to the OFF position.
2. Connect an aseptic battery housing assembled with a fully charged aseptic battery or a fully charged sterilizable battery to the handpiece.
3. The power will be engaged when the battery is installed. This is evident by the gas gauge illuminating on the aseptic battery housing or on the sterilizable battery.
4. Pull the trigger (*UPPER trigger for double trigger handpiece*) and verify the handpiece does not turn ON (no motor motion) with the MODE SWITCH in the OFF position (trigger will be blocked).
5. Pull the LOWER trigger (*double trigger handpiece only*) and verify the handpiece does not turn ON (no motor motion) with the MODE SWITCH in the OFF position (trigger will be blocked).

6. Turn the handpiece ON by sliding the MODE SWITCH to the LEFT.

7. Pull and release the trigger (*UPPER trigger for double trigger handpiece*) and verify the motor drive coupler spins and stops turning.

CAUTION: When the trigger is released, the motor will stop. If the motor continues, check the trigger to verify it returned all the way to the full out position and is not stuck partially or fully depressed. Remove from service if the motor does not stop or the trigger is not returning to the full out position.

8. Pull and release the LOWER TRIGGER (*double trigger handpiece only*) and verify the motor drive coupler spins the opposite direction as before and stops turning.

9. Fully pull both the LOWER TRIGGER and UPPER TRIGGER (*double trigger handpiece only*) and verify the motor drive coupler oscillates.

10. Change mode on the handpiece by sliding the MODE SWITCH to the RIGHT.

11. Pull and release the trigger (*UPPER trigger for double trigger handpiece*) and verify the motor drive coupler spins and stops turning.

CAUTION: The SINGLE TRIGGER HANDPIECE will enter reverse mode when the mode switch is changed. The motor will turn the opposite direction. The DOUBLE TRIGGER HANDPIECE will turn the same direction and is not reliant on the LEFT or RIGHT position of the mode switch.

12. Pull and release the LOWER TRIGGER (*double trigger handpiece only*) and verify the motor drive coupler spins (opposite direction as before) and stops turning.

13. This completes the trigger and mode switch functional verification. See the System Functions section for the full handpiece functionality.

ATTACHMENTS

ATTACHMENTS

The X Series Power System attachments have a variety of features and physical characteristics as described below.

WARNING: Before using, always verify that the attachment and cutting tool are locked into place.

ATTACHMENT INSERTION AND DISCONNECTION

1. Insert attachment by aligning teeth and pressing until “click” is heard

CAUTION: An audible click indicates that the attachment has been fully connected to the handpiece. It is recommended to test its secure fit prior to each use.

2. To test fit, lightly pull the attachment away from the handpiece to ensure that it is fully engaged
3. Disconnect by pressing / holding the button under the barrel while pulling the attachment away from the handpiece

ROTARY ATTACHMENTS

WARNING: Before using, always verify the attachment has the color coding corresponding to the desired speed, as shown in the tables below: red = 270 rpm; yellow = 600 rpm; blue = 1'500 rpm.

HIGH TORQUE REAMER ATTACHMENTS


Reference Number	Attachment Description	Rotation Speed	Illustration
89-8522-270-10	AO High Torque Reamer Attachment	270 rpm	
89-8522-270-30	Harris High Torque Reamer Attachment	270 rpm	
89-8522-270-40	Hudson High Torque Reamer Attachment	270 rpm	
89-8522-270-50	Trinkle High Torque Reamer Attachment	270 rpm	
89-8522-270-60	DCS/DHS High Torque Reamer Attachment	270 rpm	
89-8522-270-80	Zimmer/Hudson High Torque Reamer Attachment	270 rpm	
89-8522-270-90	1/4 IN High Torque Hybrid Chuck Reamer Attachment with Key (0.5 – 6.35 mm)	270 rpm	

STANDARD TORQUE REAMER ATTACHMENTS

Reference Number	Attachment Description	Rotation Speed	Illustration
89-8522-600-10	AO Standard Torque Reamer Attachment	600 rpm	
89-8522-600-20	Small AO Standard Torque Reamer Attachment	600 rpm	
89-8522-600-30	Harris Standard Torque Reamer Attachment	600 rpm	
89-8522-600-40	Hudson Standard Torque Reamer Attachment	600 rpm	
89-8522-600-50	Trinkle Standard Torque Reamer Attachment	600 rpm	
89-8522-600-60	DCS/DHS Standard Torque Reamer Attachment	600 rpm	
89-8522-600-80	Zimmer/Hudson Standard Torque Reamer Attachment	600 rpm	
89-8522-600-90	1/4 IN Standard Torque Hybrid Chuck Reamer Attachment with Key (0.5 – 6.35 mm)	600 rpm	

- All reamer attachments have a cannulation of 4 mm with the exception of the Small AO Standard Torque attachment which has a cannulation of 2.6 mm.
- For reaming procedures (270 rpm – 600 rpm), use triangular shaped shanks to avoid slippage between shank and chuck during use.
- The stated chuck diameter ranges of the hybrid chucks are specified minimums. The actual ranges of the chucks may be larger.
- The Zimmer/Hudson attachments may be referred to as “Hudson Modified Trinkle”.
- The Zimmer/Hudson attachments are compatible with Stryker Triathlon adapters and shanks.
- Triathlon is a registered trademark of Stryker Corporation. DCS/DHS are registered trademarks of Synthes Inc.

DRILLING ATTACHMENTS

Reference Number	Attachment Description	Rotation Speed	Illustration
89-8522-150-20	Small AO Drill Attachment	1'500 rpm	
89-8522-150-40	Hudson Drill Attachment	1'500 rpm	
89-8522-150-50	Trinkle Drill Attachment	1'500 rpm	
89-8522-150-80	Zimmer/Hudson Drill Attachment	1'500 rpm	
89-8522-150-90	1/4 IN Hybrid Chuck Drill Attachment with Key (0.5–6.35 mm)	1'500 rpm	
89-8522-150-91	1/4 IN Chuck Key Only		
89-8522-150-95	5/32 IN Hybrid Chuck Drill Attachment with Key (0.5–4.0 mm)	1'500 rpm	
89-8522-150-96	5/32 IN Chuck Key Only		

ATTACHMENTS

- All drilling attachments have a cannulation of 4 mm with the exception of the Small AO Drill attachment, which has a cannulation of 2.6 mm; and the 5/32 IN drilling attachment, which has a cannulation of 3 mm
- The stated chuck diameter ranges of the hybrid chucks are specified minimums. The actual ranges of the chucks may be larger
- The Zimmer/Hudson attachments may be referred to as “Hudson Modified Trinkle”
- The Zimmer/Hudson attachment is compatible with Stryker Triathlon® adapters and shanks.
- Triathlon® is a registered trademark of Stryker Corporation.

Connecting the reamer and drill attachments to the X Series Power System

1. Place the handpiece safety switch in the OFF position.

2. Hold the handpiece so that the reaming attachment is facing away from the body.

3. Slide the collet sleeve backward toward the handpiece; insert the accessory shank into the collet.

4. Release the collet sleeve to lock the accessory.

5. Pull the accessory shank to ensure that it has locked into the reaming attachment.

6. To remove the accessory, slide the collet sleeve backward toward the handpiece and pull accessory shank out of the attachment.

Connecting the Shank on the X Series Power System Hybrid Chuck attachments:

1. Place the handpiece safety switch in the OFF position.
2. Open the key/keyless attachment by rotating the collet counter clockwise by hand. Do not over open the jaws of the attachment or use the power of the handpiece to open the collet as this may jam the attachment in the open position.
3. Insert the shank into the collet and close the jaws of the attachment by rotating the collet clockwise. Ensure that the shank is centered in the jaws.
4. Tighten the collet initially by hand and then use the key if further tightening is required.
5. Tug on the shank to ensure that it has locked into the attachment.

ATTACHMENTS

6. To remove the drill, rotate the collet counter clockwise, by hand or using the key if required, and slide the drill shank out of the attachment.

SAW ATTACHMENTS

NOTE: The X Series Power System saw attachments are designed to be used with compatible saw blades sold by Zimmer Biomet.

CAUTION: Do not open the collet and/or manipulate the saw blade directly over the wound. Manipulate or change the saw blade on the back table with the handpiece turned off.

SAW ATTACHMENTS

Reference Number	Attachment Description	Speed	Illustration
89-8522-450-60	Sagittal Saw Attachment	12'500 cpm	
89-8522-450-20	Reciprocating Saw Attachment	12'500 cpm	

NOTES:

- The X Series Power System saw attachments can be connected in eight (8) different positions.
- Some X Series Power System sagittal saw blades can be attached in 5 different lengths positions. Align the arrows on the sagittal saw attachment hub with the desired blade length marking on the saw blade.

Connecting the saw blade to the X Series Power System sagittal saw attachment:

CAUTION: Follow the below instructions to ensure the saw blade is properly affixed and its length properly adjusted in the sagittal saw attachment. Handle the saw blades with great care when a power source is connected to the handpiece.

1. Rotate the knob to the counter clockwise to open/unlock.

2. Insert the saw blade into the collet to the desired position by aligning the “arrows” on the sagittal saw attachment hub with the desired length; ensure proper fit of the blade over the pins.

3. Grab the knob firmly, avoid turning with fingertips only. Rotate the knob clockwise to close/lock. Ensure that the saw blade is tightly locked by turning until the clicking sounds stop.

4. Tug on the saw blade to ensure that it is tight in the attachment.

5. To adjust the saw blade length
 - Rotate the knob counter clockwise to release the saw blade. It may be necessary to press down on the front portion of the saw blade to loosen

- Adjust the saw blade length by sliding the saw blade forward/backward to the desired position. Ensure desired length by aligning the “arrows” on the sagittal saw attachment hub with the desired length; ensure proper fit of the saw blade over the pins

- Rotate the knob clockwise to close/lock. Ensure that the saw blade is tightly locked by turning until the clicking sound stops

6. To remove the saw blade, rotate the knob counter clockwise and remove the saw blade. It may be necessary to press down on the front portion of the saw blade to loosen and then pull free.

Connecting the saw blade to the X Series Power System reciprocating saw attachment:

CAUTION: Follow the below instructions to ensure the saw blade is properly affixed in the reciprocating saw attachment. Handle the saw blades with great care when a power source is connected to the handpiece.

1. Rotate the collet counter clockwise and maintain that position.
2. Introduce the saw blade into the guide while holding the collet open. Release the collet to lock the saw blade.
3. Tug on the saw blade to ensure that the saw blade is locked in the attachment.
4. To remove the saw blade, rotate the collet counter clockwise, maintain that position, and remove the blade.

SPECIALTY ATTACHMENTS

STERNUM SAW GUARDS

The X Series Power System sternum saw guards are designed to be used with the X Series Power System reciprocating saw attachment and compatible saw blades sold by Zimmer Biomet®.

Sternum Saw Guards

Reference Number	Attachment Description	Cutting Length	Illustration
89-8522-440-90	Sternum Saw Guard – Standard	24 mm	
89-8522-440-95	Sternum Saw Guard – Large	38 mm	

NOTE: The X Series Power System sternum saw guards can be connected in two (2) different positions allowing the clinician to hold the handpiece in the direction preferred.

Connecting the sternum saw guard and saw blade to the X Series Power System reciprocating saw attachment:

CAUTION: Follow the below instructions to ensure the sternum saw guard is properly affixed to the reciprocating saw attachment.

1. Rotate the collet on the reciprocating saw attachment and maintain that position.

2. Introduce the saw blade into the attachment while holding the collet. Release the collet to lock the saw blade.

3. Lightly pull / tug on the saw blade to ensure that the saw blade is locked in the attachment.

4. In one hand, pull back the collet of the sternum saw guard as shown.

5. Slide the sternum saw guard onto the reciprocating saw attachment. Ensure that the grooves on both parts align and that the saw blade is inserted into the slot at the tip. The collet will click into place when everything is well aligned.

6. To remove the sternum saw guard, pull back the collet and slide the sternum saw guard off of the reciprocating saw attachment. If the sternum saw guard is pulled by the collet, it will automatically slide off of the reciprocating saw attachment.

7. To remove the saw blade, rotate the collet, maintain that position, and remove the blade. The saw blade can only be removed when the sternum saw guard has been removed.

PIN AND WIRE DRIVERS

The X Series Power System pin and wire driver attachments allow the surgeon to drill pin diameters of 2.0 to 4.0 mm and wire diameters of 0.6 to 2.2 mm through the handpiece.

Pin and Wire Driver Attachments

Reference Number	Attachment Description	Rotation Speed	Illustration
89-8522-451-42	Pin Driver 2.0–4.0 mm	1'500 rpm	
89-8522-451-43	Wire Driver 0.6–2.2 mm	1'500 rpm	

Procedure to use the X Series Power System pin and wire driver attachments:

1. Select the appropriate pin/wire driver attachment according to the pin/wire diameter and insert the attachment into the handpiece.

2. Insert the pin/wire in the X Series Power System pin/wire driver attachment and leave the desired length exposed.

3. Ensure the pin/wire is secured by pulling slightly on the tension lever and checking if the pin/wire is tightly locked.

4. For drilling, pull slightly on the tension lever so that the pin/wire is secure.

5. To advance the pin/wire, release the tension lever and move the X Series Power System pin/wire driver attachment in the desired direction, then pull the tension lever again. Repeat this procedure as needed.

Attachment Disposal: The X Series attachments are made with metals and polymers and should be recycled when they reach their expected service life. The attachments should be properly disposed of by a qualified recycler in order to protect the environment. Do not mix them with batteries and/or electrical wastes. Contact your local Zimmer Biomet representative or your local recycler for recycling or disposal information.

POWER SOURCES AND CHARGERS

BATTERIES

The X Series Power System batteries have a variety of features and physical characteristics as described below.

STERILIZABLE BATTERY

REF 89-8521-470-20

ASEPTIC BATTERY

REF 89-8521-470-30

Ref.	Feature	Action								
1	Battery release latch	Battery release latch that securely latches the aseptic battery housing or the sterilizable battery to the handpiece. When fully depressed, the aseptic battery housing or the sterilizable battery can be removed from the handpiece.								
2	Battery Status Lights*	Status indicators for the battery: <table><tr><td>Battery capacity at less than 25%</td><td> </td></tr><tr><td>Battery capacity at 25% – 49%</td><td> </td></tr><tr><td>Battery capacity at 50% – 74%</td><td> </td></tr><tr><td>Battery capacity at 75% – 100%</td><td> </td></tr></table>	Battery capacity at less than 25%		Battery capacity at 25% – 49%		Battery capacity at 50% – 74%		Battery capacity at 75% – 100%	
Battery capacity at less than 25%										
Battery capacity at 25% – 49%										
Battery capacity at 50% – 74%										
Battery capacity at 75% – 100%										
3	Aseptic Battery Housing or Sterilizable Battery Terminals	Aseptic Battery Housing or Sterilizable Battery terminals for transmitting energy and signals between the battery and the handpiece.								
4	Aseptic Battery Housing REF 89-8521-470-40	Sterilizable housing that is used to hold the red aseptic battery. When combined with the aseptic battery it creates the aseptic battery system used to power the handpiece.								
5	Aseptic Battery Housing Lid Release Mechanism	Release mechanism that allows the lid of the aseptic battery housing (to hold the aseptic battery) to be opened.								
6	Aseptic Battery Housing Lid	Closure mechanism for the aseptic battery								
7	Aseptic Battery REF 89-8521-470-30	Red Aseptic Battery that is inserted into the Aseptic Battery Housing via an aseptic transfer technique to create the aseptic battery system.								

Ref.	Feature	Action
8	Aseptic Battery Pull Handle	Pull handle on the aseptic battery used for pulling the battery out of the aseptic battery housing.
9	Aseptic Battery Funnel REF 89-8521-470-50	Sterilizable funnel used to transfer the aseptic battery into the aseptic battery housing
10	Sterilizable Battery REF 89-8521-470-20	Black Sterilizable Battery that connects directly to the handpiece.

*For the aseptic battery, the left picture shows the aseptic battery housing into which the aseptic battery is inserted in order to connect to the handpiece.

SETUP

During shipping and storage, the X Series Power System batteries could become weak. Primary and replacement batteries are not shipped with a full charge. Always charge the battery approximately 2 hours before initial setup, including initial checks, tests and any institution-performed biomedical evaluations.

BATTERY CHARGERS

The battery chargers have the ability to charge up to six batteries at once and contain the features and physical characteristics described below.

CHARGER FOR ASEPTIC BATTERY

REF 89-8521-470-00

CHARGER FOR STERILIZABLE BATTERY

REF 89-8521-470-10

Ref.	Feature	Action
1	Battery Charger	Six bay charger used to charge the aseptic battery or the sterilizable battery.
2	Charging Dock Display	This display is used to provide information to the user relative to the battery charger function and the battery charging status. There is an indicator on each bay of the battery charger that operates independent of each of the others.
3	Blue Status Indicator	This indicator is used to provide information to the user relative to charging circuits and battery conditions. This is accomplished through an illumination or flashing. There is an indicator on each bay of the battery charger that operates independent of each of the others.
4	Rear power entry module switch and fuse holder	Power entry module switch on the back of the battery charger used to set the charger to OFF. The power entry module contains fuses for the AC Mains.
5	Potential Equalization Conductor Stud	Conductor stud that can be used to provide a connection between the battery charger and the potential equalization bus-bar of the electrical installation.

POWER SOURCES AND CHARGERS

SETUP

The battery charger for aseptic battery will charge the aseptic battery only.

The battery charger for sterilizable battery will charge the sterilizable battery only.

WARNING: Risk of electric shock. Do not remove cover.

WARNING: Never use any charger other than an X Series Power System battery charger to charge the X Series Power System battery.

CAUTION: To ensure the battery is properly charged, insert the battery and check the indicator lights as described below.

NOTE: The pictures below are showing the charger for aseptic battery (red battery) which is similar to the charger for sterilizable battery (black battery).

1. Attach the power cord at the rear of the battery charger in the power entry module.

2. The battery charger can be powered ON before or after the batteries are inserted.

3. The blue status indicator will flash when the battery charger is starting. The charging dock display will show the Startup Sequence pictogram when it is starting (outline of a battery cell with a clock inside). Once the startup sequence is complete, the status indicator will be ON, the charging dock display will be not illuminated and ready to charge when no battery is present.

4. If not inserted, slide a battery into the battery charger until the battery is fully engaged. The battery terminals are face down and the battery status LEDs are inserted into the battery charger.

5. The blue status indicator will flash when the battery is initially inserted to indicate communication is being established. The LED status indicator will stay ON while the charging dock display will flash with a battery cell outline to indicate the communication is establishing.

6. The blue status indicator will be ON when the battery communication has been established and the battery is charging. The charging dock display will show the charging pictogram: the outline of a battery cell and four charging bars. The charging bars will flash according to the charging status of the battery.

7. The charging dock display will show the status of the battery charge capacity. This display will also show fault status if a fault is found.

8. During the charging, the charging dock display will show the battery charge capacity status. The blue status indicator will be flashing during this period.

Charge Icon	Action	Blue LED	Battery Charge Level
(Blank)	OFF	ON	N/A
	First bar flashing.	ON	0% to 24%
	First bar continuously displayed and second bar flashing.	ON	25% to 49%
	Two first bars continuously displayed and third bar flashing.	ON	50% to 74%
	Three first bars continuously displayed and fourth bar flashing.	ON	75% to 95%
	Steady	ON	96% to 100%

9. If the battery charger circuit encounters a fault, the charging dock display will display the error code. The following error codes could be shown depending on the fault:

Charge Icon	Action	Blue LED	Root Cause	Fault Identifier
	Alternating	Flashing	Input voltage fault	10
	Alternating	Flashing	Output voltage fault	20
	Alternating	Flashing	PCB over temperature fault	30
	Alternating	Flashing	PCB under temperature fault	40
	Alternating	Flashing	Cooling fan fault	50
	Alternating	Flashing	Internal clock fault	60
	Alternating	Flashing	Output current fault	80
	None	Flashing	Communication fault	None

POWER SOURCES AND CHARGERS

10. If the battery charger encounters a battery fault, the charging dock display will display the error code. The following error codes could be shown depending on the fault:

Charge Icon	Action	Blue LED	Root Cause	Fault Identifier
	Alternating	Flashing	Communication fault	05
	Alternating	Flashing	Authentication fault	15
	Alternating	Flashing	Over current fault	25
	Alternating	Flashing	Over voltage fault	35
	Alternating	Flashing	Battery over temperature fault	45
	Alternating	Flashing	Battery under temperature fault	55
	Alternating	Flashing	Over charge current fault	65
	Alternating	Flashing	Charge timeout fault	75
	Alternating	Flashing	Open battery fault	85
	Alternating	Flashing	Pre-charge time fault	95

11. Once the battery has been successfully charged, allow to cool to room temperature before using.

WARNING: Never attempt to charge a battery that has obvious signs of damage. Remove a battery from service if there are any signs of damage or sterilization abuse.

CAUTION: Never charge or discharge a battery that is hot. The battery charger and/or handpiece will test the battery temperature and will not allow a hot battery to be charged or discharged.

CAUTION: Sterilizing a hot battery may cause the battery cells to be exposed to higher internal temperatures. Higher internal temperatures may shorten the battery's overall expected service life.

FUNCTIONALITY CHECK

This section describes automatic and manual checks to verify system functions. The diagnostics are built into the battery charger and batteries and will automatically check the function upon power-up. The manual checks can be performed by following the instructions below.

1. Verify the battery charger is plugged into ~ AC power mains.
2. Turn the power ON by toggling the power switch at the rear of the battery charger to ON.
3. The battery charger will execute a startup sequence. The battery charger will show a battery clock icon indicating the startup sequence. The blue LED will be flashing.

4. If there is no battery present in the bay, the display will be OFF. The display will automatically illuminate when a battery is inserted. The blue LED will remain ON.
5. Next, insert a battery into the battery charger. The display will flash an empty battery icon indicating it is communicating with the battery and reading status. The charging display will be flashing.

6. The battery charger will analyze the battery automatically for faults. If no faults are detected, the battery will begin to charge. The charging icon being displayed depends upon the capacity of charge already in the battery.

7. If a fault is detected within the battery charger, the blue LED will be flashing. The icon will alternate between the fault icon and a fault code. A Fault 10 is shown in the example below.

8. If a fault is detected within the battery, the blue LED will be flashing. The icon will alternate between the fault icon and the fault code inside a battery outline. A Fault 05 is shown in the example below.

9. Inspect all bays and ensure they are working as described. Inspect all batteries and ensure they are working as described.
10. Attach the battery to a known functional handpiece.
11. Verify the blue battery status lights illuminate on the battery. Pull the handpiece trigger and verify the motor turns ON.
12. Disconnect the battery.
13. This completes the battery charger and battery functional verification. See the System Functions section for the full battery charger and battery functionality.

WARNING: The batteries cannot be returned by air transport (in accordance with IATA DGR). Only GROUND shipment is allowed for these products. For additional information, please contact your local Zimmer Biomet representative.


Li-ion

Battery Disposal: The battery is of a Lithium-Ion type. When the battery reaches the end of its useful life, the battery should be disposed of by a qualified recycler or hazardous materials handler. Do not mix this battery with the solid waste stream. Contact your local recycler or hazardous material handler for recycling or disposal information.

STERILIZATION BASKETS

STERILIZATION BASKET FOR USE WITH ASEPTIC BATTERY HOUSINGS

- | | | |
|------------------------|---------------------------------|----------------------------|
| 1. Handpiece (2 Slots) | 4. Rotary Attachments | 7. Sagittal Saw Attachment |
| 2. Pin Driver | 5. Reciprocating Saw Attachment | 8. Aseptic Battery Funnel |
| 3. Wire Driver | 6. Hybrid Chuck Key | 9. Aseptic Battery Housing |

SYSTEM FUNCTIONS

OPERATION

WARNING: To isolate the device from ~ AC Power Mains in case of an emergency, disconnect the power cord from the wall outlet or from the appliance inlet at the rear of the device. Do not position the device in such a way that obstructs the access to the wall outlet or appliance inlet.

SURGEON DECISIONS

The surgeon's discretion will be used to determine the following:

- When to use the X Series Power System.
- Desired power source.
- What attachments are to be used.
- The speed that is appropriate for the different rotating cutting tools.
- The shape, length, width, and diameter of the cutting tools.

The X Series Power System has three primary modes: Forward, Reverse and Oscillate.

- Forward – Drives the reaming and drilling attachments in the forward cutting direction. This also drives the sagittal and reciprocating saw attachments in the cutting direction.
- Reverse – Drives the reaming and drilling attachments in the reverse non-cutting direction. This also drives the sagittal and reciprocating saw attachments in the cutting direction.
- Oscillate – Drives the reaming and drilling attachments in the forward cutting direction and then in the reverse non-cutting direction. This function is only available on the double trigger handpiece.

CAUTION: The oscillating mode is only intended for drilling function. Oscillate mode is not intended to be used for the sagittal or reciprocating saw attachments. The performance of the saw attachments and the battery life will be reduced if attempting to use a saw function in oscillate mode.

TESTS AND CHECKS

The X Series Power System components shall produce the results explained in the following visual inspection steps and in the aforementioned functional testing steps exactly as indicated. Failure to do so indicates that a problem may exist and the device is not to be used until necessary repair has been made. It is strongly recommended to refer to Maintenance section.

It is strongly recommended that the visual inspection and functional checks of the handpiece and battery charger be performed every three (3) months or after any scheduled or unscheduled maintenance. Failure to perform routine visual inspection and functionality checks may result in undesired performance.

CAUTION: If the X Series Power System fails to produce the expected results as described then the affected device must be removed for service.

MAINTENANCE INDICATORS

- This section is intended to assist the user in determining whether a X Series Power System instrument has worn to an extent that it is no longer suitable for use. For additional reusable instrument care instructions, see the respective specific sections of the Owner's Manual.
- Several indicator categories that could render the instrument no longer suitable for use are described below. These indicators are signs of wear and damage.
- If any of these indicators is observed on an instrument, contact Zimmer Biomet customer service or your local representative and return the product for investigation and repair if applicable.
- Not following this instruction can lead to an unusable device during surgery that could potentially lead to harm to user or patient. All visual assessment of indicators shall be executed after cleaning and sterilization of the devices.

INDICATOR CATEGORY	DESCRIPTION	SPECIFIC CHECKS
Fracture	Cracked or broken into two or more pieces.	Check especially battery housing and all other polymer components.
Deformation	Bent, deformed, melted.	Check especially the battery housing and perform visual check of seal(s). Check especially the distal part of the sternum guard for deformation.
Surface Damage	Wear, excessive scratches/dents, flaking, abrasion of top layer or coating.	Check especially crown body and other mating surfaces.

SYSTEM FUNCTIONS

INDICATOR CATEGORY	DESCRIPTION	SPECIFIC CHECKS
Discoloration	Loss of color, functional color codings (attachment color rings).	Check especially surface of the device and color coding of the attachments.
Corrosion	Potential corrosion or traces of corrosion on the device.	Check especially in openings or undercuts and grooves.
Marking Damage	Marking and/or UDI unreadable, unrecognizable or damaged (including barcodes).	Check markings on the device.
Leakage	Liquid or other substances leaking from device. Cracked or damaged seals, without visible leakage shall also trigger a return of the device.	Check especially rotating parts and seal surfaces.
Functional damage	Does not work as intended (e.g does not turn, lock, open, close, assemble, disassemble or otherwise function as intended, loosening, vibration, absence of oscillating mode, incorrect directionality of motor, missing or degraded speed ramp, loss of speed control) Non-function of the device.	Refer to the individual component sections in the Owner's Manual for specific functional check.
Particle Release	Debris, particle release.	Check especially in openings or undercuts and grooves. Check especially for particles falling.
Overheating	Device too hot to handle during operation.	If users experience pain or discomfort due to elevated temperatures of the device, the device shall be returned.

CAUTION: Prior to initial usage, Zimmer Biomet X Series Power System must be cleaned and sterilized per the instructions (See Cleaning, Decontamination and Sterilization Instructions).

WARNING: Before attaching a power source to the handpiece, always slide the handpiece mode switch to the center OFF position.

ASSEMBLY OF THE HANDPIECE, ATTACHMENT AND BATTERIES

Both the double trigger and single trigger handpieces will accept all X Series Power System attachments and power sources. The attachment will be application-specific.

1. Before attaching a power source, verify the mode switch is set to the center position. This will lock the trigger to avoid unintended activation.
2. Follow aseptic transfer protocol when passing sterile materials between the scrubbed person and the circulating person.

WARNING: The aseptic transfer protocol must be followed without any deviation in order to maintain the sterility of the aseptic battery housing. The aseptic transfer must only be executed using a sterile aseptic battery housing and aseptic battery funnel.

ASEPTIC BATTERY TRANSFER

- a) **Scrubbed person:** open the aseptic battery housing lid by pressing on both metallic latches simultaneously.

- b) **Scrubbed person:** assemble the aseptic battery funnel with the aseptic battery housing.

NOTE: When seated correctly, the aseptic battery funnel will cover the entire edge of the aseptic battery housing opening.

- c) **Circulating person:** insert the aseptic battery into the aseptic battery housing through the aseptic battery funnel. The pull handle can be utilized to hold the battery (and to remove it during disassembly prior to the post-operative reprocessing). Ensure the battery is fully seated in the aseptic battery housing and then remove the aseptic battery funnel.

WARNING: Avoid all contact with the aseptic battery housing. If the circulating person or the aseptic battery contacts the outside of the aseptic battery housing, it must be cleaned and resterilized before any further use in the operating room.

- d) **Scrubbed person:** push the aseptic battery housing lid shut by contacting only its outside face. Lock the lid by pushing it until both metallic latches snap in place.

WARNING: A double click indicates that both metallic latches are properly closed. This should also be verified visually.

CLEANING, DECONTAMINATION AND STERILIZATION

WARNING: Avoid all contact with the inside of the aseptic battery housing and with the aseptic battery. If the scrubbed person contacts the aseptic battery, he/she must rescrub.

e) The **scrubbed person** can now attach the aseptic battery housing on to the handpiece.

CLEANING, DECONTAMINATION AND STERILIZATION

These instructions cover the entire Zimmer Biomet X Series Power System and its accessories.

CAUTION: The battery chargers and aseptic battery are NOT washable or sterilizable. They are designed to be cleaned with a cloth that has been dampened with a mild detergent.

CAUTION: The X Series Power System components are not intended to be washed using detergents with pH exceeding 10.8.

CAUTION: For automatic cleaning and sterilization, the reciprocating saw attachment shall be oriented such that the tip points downward and the handpiece connector points upwards. As a secondary option, it can be positioned horizontally.

CAUTION: The Sterilizable battery is not intended to be cleaned in an automatic washer.

WARNING: The X Series Power System components including the batteries and battery chargers are not intended to be immersed in a cleaning solution or cleaned in an ultrasonic bath.

WARNING: Use of the Zimmer Biomet sterilization case is strongly recommended. The processes described in this manual are the only cleaning, disinfection and sterilization processes covered by the warranty. Processes outside of the parameters may decrease performance and/or reliability.

WARNING: Final responsibility for verifying cleaning and sterility, using the equipment and processes of the healthcare facility, and the equipment and parameters supplied by Zimmer Biomet, lies with the healthcare facility. To ensure optimal processing all cycles and methods should be validated for different sterilization chambers, wrapping methods and/or various load conditions.

WARNING: When placing components into a basket for cleaning or sterilization, avoid all contact between the devices. Touching devices may cause damage during washing and may jeopardize the efficacy of the reprocessing. Ensure cannulations point downwards to facilitate drying of the instruments following automatic cleaning.

WARNING: All system components and cutting tools must be disconnected from each other prior to reprocessing in order to ensure the efficacy of the cleaning, decontamination and sterilization processes.

These instructions must be performed by professional staff appropriately trained in cleaning and disinfection methods.

PRECAUTIONS FOR CLEANING, DECONTAMINATION AND STERILIZATION

- Always follow universal precautions for infection control and treat all human blood and certain human body fluids as if they were known to be infectious.
- It is highly recommended that the X Series Power System be cleaned within 30 minutes of use to minimize the chance for biological fluids and materials to dry.
- Consult the cleaning solution manufacturer's instructions prior to use to ensure the cleaning solutions are compatible with aluminum and stainless steel.
- Do not use cleaning solutions with PHENOL, ALDEHYDE or CHLORINE BLEACH.
- Do not immerse handpiece or batteries in any liquids.
- Remove the attachment from the handpiece prior to cleaning.
- Remove the sternum guard from the reciprocating saw attachment prior to cleaning.
- Remove the sterilizable battery or the aseptic battery housing from the handpiece prior to cleaning.
- Not for use in ultrasonic cleaning systems.
- Repeat washing at high pH can lead to restricted movement of the tool release mechanism. If the mechanism is observed to not move freely then a sterilizable instrument lubricant can be used to restore free movement. Hinge-Free Instrument Lubricant from STERIS has been shown to be effective at restoring movement. Lubricant intended for steel medical instruments that will not impeded the sterilization method used at your facility can be used. Always follow your hospital protocol for automated or manual lubrication to restore free movement. Hinge-Free is a registered trademark of STERIS.

CAUTION: To prevent corrosion, do not use a cleaning solution with a pH over 10.8.

PREPARING FOR MANUAL CLEANING

Immediately after use, carefully rinse the device with reverse osmosis (RO) or distilled water in order to remove all visible soil. If the cleaning cannot be done immediately, put the disconnected devices in a new soft cloth heavily saturated with cleaning solution per hospital protocol to prevent drying.

1. Place the MODE SWITCH in the OFF position.
2. Disconnect the sterilizable battery or the aseptic battery housing from the handpiece by pressing the blue tab and gently pulling the sterilizable battery or the aseptic battery housing toward the rear (away from handpiece).

CAUTION: The saw blades and drill bits are sharp. Avoid any contact with the cutting tools when a power source (battery) is attached to the handpiece.

3. Disconnect the single use saw blade or drill bit and dispose of properly.
4. Disconnect the attachment from the handpiece.

CAUTION: If a sternum saw guard was used then disconnect it from the reciprocating saw attachment by first pulling back the collet sleeve on the sternum saw guard.

5. If using the aseptic battery, open the lid of the aseptic battery housing by pressing on both metallic latches simultaneously and remove the battery. If needed, use the metallic handle on the battery to pull it out of the aseptic battery housing.

MANUAL CLEANING MATERIALS

1. Cold/room temperature RO water or distilled water with a hardness lower than 120 mg/L of CaCO₃.
2. Detergent with pH not exceeding 10.8, to be diluted using the RO or distilled water according to the cleaner manufacturer's recommendations.
3. Non-abrasive brush.
4. Non-abrasive cannulation brush for cleaning handpiece and attachment cannulation.
5. Non-shedding dry cloth.

CLEANING, DECONTAMINATION AND STERILIZATION

MANUAL CLEANING PROCESS

Manual cleaning is required to remove gross debris and traces of blood, debris and stains from usage. This is a critical step in the cleaning, decontamination and sterilization strategy.

Step	Phase	Directions
Step 1	<i>RINSE</i>	Rinse the handpiece, attachments, sterilizable battery and aseptic battery housing with RO or distilled water until all visible body fluids and tissues have been removed.
Step 2	<i>SCRUB</i>	For the handpiece: Wrap the handpiece in a cloth that has been soaked in detergent and wait for 10 minutes to let the detergent work. Alternatively, spray the handpiece with an enzymatic cleaning foam and wait for the amount of time indicated by the foam supplier in order to let the foam work. Then, carefully clean each part of the handpiece, attachments, sterilizable battery and aseptic battery housing with a soft brush dampened with the appropriate detergent. Scrub for a minimum of 2 minutes paying particular attention to critical areas such as exposed parts, narrow cavities, small openings and moving parts. Let the detergent flow into the cavities and areas that may be more difficult to reach. Use an instrument brush to clean the cannulations of the handpiece and attachments.
Step 3	<i>RINSE</i>	Rinse carefully with RO or distilled water a minimum of 30 seconds so that all the detergent is fully removed.
Step 4	<i>DRY</i>	Dry the handpiece, attachments, sterilizable battery and aseptic battery housing using a non-shedding cloth or medical quality filtered compressed air.

WARNING: Not following the above manual cleaning steps may result in insufficient removal of biological contaminants or of cleaning agents.

MANUAL CLEANING INSPECTION

Inspection after the manual cleaning is important to ensure all visible soil and detergent is fully removed. If any visible soil is noticed, repeat the RINSE, SCRUB, RINSE, DRY process.

WARNING: If the devices are not disinfected immediately after rinsing, carefully dry the device with a non-shedding cloth to avoid any microbial contamination due to humidity.

DISINFECTION

The disinfection process provides a method to disinfect the handpiece, attachments, sterilizable battery and aseptic battery housing after inspection has shown the device is clean.

Step	Phase	Directions
Step 1	<i>COVER</i>	Cover the device in a soft cloth heavily saturated with the recommended detergent/disinfectant for 30 minutes.
Step 2	<i>WIPE</i>	After the 30 minute contact wait time, wipe the device down with the saturated cloth, paying particular attention to joints, mating areas, crevices, and small openings.
Step 3	<i>RINSE</i>	Rinse carefully with RO or distilled water so that all the detergent is fully removed.
Step 4	<i>DRY</i>	Dry the handpiece, attachments, sterilizable battery and aseptic battery housing using a non-shedding cloth or medical quality filtered compressed air.

PREPARING FOR AUTOMATED CLEANING

CAUTION: The Sterilizable battery is not intended to be cleaned in an automatic washer.

Automated cleaning following manual cleaning is recommended but not mandatory.

After performing the Manual Cleaning, Inspection and Disinfection process above, inspect the handpiece, attachments and aseptic battery housing once again for any signs of residual soil. Repeat the entire Manual Cleaning, Inspection and Disinfection process if signs of residual soil are detected.

AUTOMATED CLEANING MATERIALS

1. Detergent with pH not exceeding 10.8.
2. Rinse aid.
3. Automated washing / sanitizer system.

AUTOMATED CLEANING PROCESS

The Zimmer Biomet basket has been validated and is recommended for the automated cleaning and disinfection process. Place all devices into their assigned positions. Generally, instruments shall be placed into the basket in such a way that cannulations point downwards, the instruments do not contact each other and the basket is not overloaded. This helps ensure the instruments are cleaned and dried properly. Then put the entire basket into the automated washer.

Establish the following process parameters in the automated washer:

Step	Phase	Temperature	Time
Step 1	PRE-WASH	20 °C (68 °F)	2:00 minutes
Step 2	CLEANING	55 °C (131 °F)	5:00 minutes
Step 3	NEUTRALIZING	20 °C (68 °F)	2:00 minutes
Step 4	RINSE I	20 °C (68 °F)	2:00 minutes
Step 5	FINAL RINSE	93 °C (200 °F)	5:00 minutes
Step 6	DRY	99 °C (210 °F)	35:00 minutes

WARNING: Only the washing process parameters shown above have been validated for cleaning effectiveness. The performance and/or reliability are not affected by lower temperatures and lower dwell times. However, the performance and/or reliability may decrease with increased temperatures and increased dwell times.

AUTOMATED CLEANING INSPECTION

Inspection after the automated cleaning is important to ensure all visible soil and detergent is fully removed. If any visible soil is noticed, repeat the PRE-WASH, CLEANING, NEUTRALIZING, RINSE I, FINAL RINSE, DRY process after the device has cooled to room temperature.

PREPARING FOR STERILIZATION

After performing the prescribed Cleaning, Inspection and Disinfection processes above, inspect the handpiece, attachments, sterilizable battery and aseptic battery housing once again for any signs of residual soil. Repeat the cleaning process if signs of residual soil are detected.

The referenced sterilization cycle parameters listed in this manual have been validated and provide a sterile result by steam sterilization according to ANSI/AAMI/ISO 17665-1:2006, Sterilization of Health Care Products. Part 1: Moist Heat. The steam sterilization process for these products has been validated to provide a sterility assurance level (SAL) of 10^{-6} or better when using the prescribed sterilization case.

Packaging used in the sterilization of these devices shall: follow EN ISO 11607; be suitable for the chosen sterilization method; provide sufficient protection of the device as well as protection of the packaging to prevent mechanical damage; be the appropriate type/grade for weight of device.

Charging of the sterilizable battery must occur prior to the sterilization step. Refer to the battery chargers section of this manual for the charging procedure. The exception to this is when performing Immediate Use Steam Sterilization without recharging the battery.

WARNING: There are different validated sterilization processes for the X Series Power System components. Follow the custom sterilization cycle parameters defined in the following sections for the different components.

WARNING: Do not use hot air, ethylene oxide (EtO), plasma, Gamma or formaldehyde for sterilization unless explicitly permitted elsewhere in this document. These are not validated processes and can damage the equipment.

WARNING: When wrapping devices prior to sterilization, follow the relevant industry standards (e.g. ANSI/AAMI ST79) and ensure the sterilization packaging has received the necessary regulatory clearance.

STERILIZATION OF HANDPIECE, ATTACHMENTS, ASEPTIC BATTERY HOUSING AND ASEPTIC BATTERY FUNNEL

STEAM STERILIZATION

The steam sterilization procedure must be carried out in a qualified autoclave. All components should be positioned in the cases following the diagrams printed on the cases. Do not overload the cases.

The handpiece, attachments, aseptic battery housing and aseptic battery funnel are to be sterilized using a custom steam sterilization cycle defined in this manual. The manual cleaning and disinfection steps must be followed prior to sterilization regardless of the sterilization exposure parameters being used. Use only the Zimmer Biomet specified sterilization tray for processing.

1. Ensure the handpiece, attachments, aseptic battery housing and aseptic battery funnel have been cleaned and disinfected using the cleaning process.
2. Place all components in the sterilization tray per the layout guide printed in the tray.
3. Wrap per hospital protocol.
4. Sterilize using one of the validated cycles with the following process parameters:

Optional Cycle	Steam cycle type	Temperature	Steam Dwell Time	Minimum Dry Time*
1	Pre-Vacuum	132 °C (270 °F)	4 minutes	45 minutes
2	Pre-Vacuum	132 °C (270 °F)	8 minutes	45 minutes
3	Pre-Vacuum	134 °C (273 °F)	3 minutes	45 minutes
4	Pre-Vacuum	134 °C (273 °F)	18 minutes	45 minutes
5	Gravity	132 °C (270 °F)	15 minutes	45 minutes

*Dry time is highly dependent on sterilization load size. Dry time may need to be increased in 5-minutes increments for larger loads.

5. Allow to cool to a comfortable operating temperature and then remove the package from the sterilizer. Inspect for visible moisture or dampness. If the package shows signs of moisture or dampness, the tray should be repackaged and re-sterilized.
6. Because dry time can be influenced by elements such as load size, autoclave performance, wrapping, steam inflow, and cool-down function, consider increasing the dry time by 5 minutes increments to improve drying as needed.

IMMEDIATE USE STEAM STERILIZATION (IUSS)

Sterilizable devices can be sterilized using an IUSS cycle in situations intended for immediate use. This method is not intended for devices that need to be stored for future use, nor held from one case to another. The cleaning and disinfection steps must be followed prior to sterilization regardless of the sterilization exposure parameters being used.

1. Ensure the device has been cleaned and disinfected using the cleaning process.
2. Place the device in a basket. Do not wrap or pouch the device.
3. Sterilize using one of the validated cycle with the following process parameters:

Optional Cycle	Steam cycle type	Temperature	Steam Dwell Time	Minimum Dry Time
1	Pre-Vacuum	132 °C (270 °F)	4 minutes	No dry time
2	Pre-Vacuum	134 °C (273 °F)	3 minutes	No dry time

WARNING: Failure to allow the package to cool at the end of the sterilization cycle may decrease performance, compromise sterility or result in burnt patient tissue or healthcare staff.

CAUTION: Do not perform rapid cooling techniques on the load as this could damage the handpiece.

STERILIZATION OF STERILIZABLE BATTERY

STEAM STERILIZATION

The steam sterilization procedure must be carried out in a qualified autoclave.

The sterilizable battery is to be sterilized using a custom steam sterilization cycle defined in this manual. The manual cleaning and disinfection steps must be followed prior to sterilization regardless of the sterilization exposure parameters being used.

1. Ensure sterilizable battery has been cleaned and disinfected using the cleaning process.
2. Place the sterilizable battery in a basket. Alternatively, the sterilizable battery may be pouched or wrapped as per hospital protocol.
3. Sterilize using one of the validated cycle with the following process parameters:

Optional Cycle	Steam cycle type	Temperature	Steam Dwell Time	Maximum Dry Time*
1	Pre-Vacuum	132 °C (270 °F)	4 minutes	4 minutes
2	Pre-Vacuum	134 °C (273 °F)	3 minutes	4 minutes

*To ensure continued functioning of the battery, the maximum dry time must not be exceeded.

4. Allow to cool to a comfortable operating temperature and then remove the package from the sterilizer. Inspect for visible moisture or dampness. If the package shows signs of moisture or dampness, the tray should be repackaged and re-sterilized.

IMMEDIATE USE STEAM STERILIZATION (IUSS)

Sterilizable batteries can be sterilized using an IUSS cycle in situations intended for immediate use. This method is not intended for sterilizable batteries that need to be stored for future use, nor held from one case to another. The cleaning and disinfection steps must be followed prior to sterilization regardless of the sterilization exposure parameters being used.

1. Ensure the sterilizable battery has been cleaned and disinfected using the cleaning process.
2. Place the sterilizable battery in a basket. Alternatively, the sterilizable battery may be pouched or wrapped as per hospital protocol.
3. Sterilize using one of the validated cycle with the following process parameters:

Optional Cycle	Steam cycle type	Temperature	Steam Dwell Time	Minimum Dry Time
1	Pre-Vacuum	132 °C (270 °F)	4 minutes	No dry time
2	Pre-Vacuum	134 °C (273 °F)	3 minutes	No dry time

4. Allow the device to cool a minimum of 5 minutes at room temperature.

WARNING: Failure to allow the package to cool at the end of the sterilization cycle may decrease performance, compromise sterility or result in burnt patient tissue or healthcare staff.

CAUTION: Do not perform rapid cooling techniques on the load as this could damage the sterilizable battery.

CAUTION: A battery that is too hot will not be allowed to function by the handpiece until it has reached the proper internal temperature.

STERRAD

The sterilizable battery can be sterilized using STERRAD defined in this manual.

WARNING: The sterilizable battery is the ONLY component of the X Series Power System that is to be sterilized using STERRAD.

For details on the STERRAD technology and the processing guide, visit: <http://www.sterradsterilityguide.com>

STERRAD is a registered trademark of Johnson & Johnson.

V-PRO

The sterilizable battery can be sterilized using V-PRO defined in this manual.

WARNING: The sterilizable battery is the ONLY component of the X Series Power System that is to be sterilized using V-PRO.

For details on the V-PRO technology and the processing guide, visit: <https://www.steris.com/products/vpro>

V-PRO is a registered trademark of STERIS.

MAINTENANCE

GENERAL MAINTENANCE INFORMATION

While the X Series Power System has been designed and manufactured to high industry standards, it is recommended that regular cleaning, decontamination, sterilization and inspection be performed to ensure continual safe and effective operation. This section contains information to assist in the effort as well as serve as a guide to expediting unscheduled maintenance.

SERVICING

PRECAUTIONS IN SERVICE

- Rough shocks/handling should be avoided.
- Do not attempt to disassemble the equipment. It is a factory sealed unit with no user serviceable parts inside. No modification of the equipment is allowed.
- Zimmer Biomet is not liable for any device malfunction resulting from repairs or service not performed by a Zimmer Biomet authorized service center.
- **Always decontaminate the system before returning to Zimmer Biomet.**

The X Series Power System must be returned every 6 months for inspection and preventive maintenance. Annual factory calibration checks and preventive maintenance are required to verify continued accuracy.

This manual includes a blank test and evaluation report. This is provided at the end of the manual. This form may be photocopied and used to document the visual, electrical safety and functional tests.

WARNING: No modification of this equipment is allowed. The X Series Power System handpiece, batteries and battery charger contain no user-serviceable parts inside. All internal parts and components must be serviced by an authorized and qualified biomedical engineer or other person thoroughly familiar with electronic medical devices. Do not attempt to dismantle, modify or repair internal components. When testing, it is the user's responsibility to ensure that the appropriate and proper measuring equipment is used. Improper or poorly maintained equipment may result in incorrect or misleading results.

MANDATORY ACTION: If the electrical tests described in the following section cannot be performed, it is recommended that the device be returned to Zimmer Biomet for the electrical safety checks.

PROTECTIVE EARTH RESISTANCE TESTING


CAUTION: The battery charger must meet the protective earth resistance requirements before performing any additional electrical safety checks.

The protective earth resistance measurement should be conducted with the power cord attached to appliance inlet. The measuring device used must be able to deliver at least 200 mA into 500 mΩ. The open circuit voltage shall not exceed 24 V. Consult the test equipment manufacturer to ensure these specifications are met. Measure the protective earth resistance between the protective earth pin (ground pin) of the power cord and the exposed metal chassis at the back of the device (the protective equalization conductor stud may also be used). The resistance reading should not exceed 300 mΩ (milliohm).

The integrity of the earth conductor of the power cord should be performed during normal maintenance intervals. During the protective earth resistance testing, the power cord should be flexed along its entire length. If during flexing, changes in resistance are observed, the power cord's protective earth conductor is presumed to be damaged and the power cord should be immediately replaced.

EARTH (EQUIPMENT) LEAKAGE CURRENT TESTING

The earth leakage current is one recommended leakage current test for the battery charger. The maximum allowable leakage current is 10 mA in a single fault condition.


The leakage current is tested from the protective earth ground of the battery charger to exposed metal as well as the plastic enclosure.

INSULATION RESISTANCE TESTING

CAUTION: The battery charger must be disconnected from supply mains when measuring the insulation resistance.

The insulation resistance testing is conducted by measuring the Mains Part to the Protective Earth ground. This test is optional since it requires a specialized insulation resistance tester. However, the test shall be conducted if the residual current device has tripped several times, or there was liquid spilled over the device or other situations where the internal insulation properties (creepage/clearance) is suspect. The test equipment used shall be capable of performing the measurements with 500 VDC. The insulation resistance will be a very high level since the insulation is intended to provide no means for electrical conductivity. Therefore, values $\geq 70\text{M}\Omega$ (megaohm) will be considered acceptable reference value. Measurements that show a significant decrease from the previous measured value indicates a problem.

BATTERY SERVICE

The X Series Power System batteries are only intended to be used with the X Series Power System handpieces and battery chargers. The handpiece will not allow a non-Zimmer Biomet battery to be used. The battery does not contain any user serviceable components.

CAUTION: The batteries cannot be returned by air transport (in accordance with IATA DGR). Only GROUND shipment is allowed for these products. For additional information, please contact your local Zimmer Biomet representative.

BATTERY-HANDLING RECOMMENDATIONS

- Do not dismantle, open or shred the battery
- Do not expose the battery to heat or fire. Avoid storage in direct sunlight.
- Do not short-circuit the battery. Do not store batteries haphazardly in a box or drawer where they may short-circuit each other or be short-circuited by other metal objects.
- Do not remove a battery from its original packaging until required for use.
- Do not subject a battery to mechanical shock.
- In the event of a cell leakage, do not allow the liquid to come in contact with the skin or eyes. If contact has been made, wash the affected area with copious amounts of water and seek medical advice.
- Do not use any charger other than that specifically provided for use with the equipment.
- Do not use any battery which is not designed for use with the equipment.
- Seek medical advice immediately if a cell or a battery has been swallowed.
- Always purchase the correct battery for the equipment.
- Keep batteries clean and dry.
- Secondary batteries need to be charged before use.
- After extended periods of storage, it may become necessary to charge and discharge the battery several times to obtain maximum performance.
- Secondary batteries will give their best performance when they are operated at normal room temperature ($68^{\circ}\text{F} \pm 9^{\circ}\text{F}$ / $20^{\circ}\text{C} \pm 5^{\circ}\text{C}$).
- Retain the original product literature for future reference.
- Use only the battery in the application for which it was intended.
- When possible, remove the battery from the equipment when not in use.
- Dispose of properly. See battery disposal instructions for further guidance on the batteries.

Battery Disposal: The battery is of a Lithium-Ion type. When the battery reaches the end of its useful life, the battery should be disposed of by a qualified recycler or hazardous materials handler. Do not mix this battery with the solid waste stream. Contact your local recycler or hazardous material handler for recycling or disposal information.

FUSE AND FUSE DRAWER REPLACEMENT

CAUTION: Use only the size and rating of the fuses described in the CHARACTERISTICS section.

To remove and install a fuse drawer in the battery charger:

- Locate the power entry module at the rear of the battery charger where the power cord attaches.
- If still attached, remove the power cord to ensure no power is being supplied to the battery charger.
- Using a tool such as a slotted small screwdriver, gently pry the holding clip of the fuse drawer to release the latch. The fuse drawer and fuses will release and can be removed.
- Push the fuse drawer into the power entry module until the fuse drawer locks into place (an audible click will be heard).
- Reconnect the power cord to power entry module, plug in to AC and verify battery charger AC power is active.

UNSCHEDULED MAINTENANCE

The X Series Power System is designed with several specific self-test features to assist in fault isolation. In the unlikely event the system malfunctions, the device is designed to stop in a SAFE STATE. This SAFE STATE will keep the motor from turning. The system should be removed from service immediately along with all accessories and reported to the designated individual responsible for equipment repair within the healthcare organization. The system should be returned to Zimmer Biomet for repair.

TROUBLESHOOTING GUIDE

To aid in unscheduled maintenance, the **TROUBLESHOOTING** table delineates a number of possible malfunctions that could occur with the unit. The most likely causes are shown for each symptom. While it is not practical to enumerate every conceivable malfunction and all possible causes, the table will assist in isolating the most common problems.

VOLTAGE ELECTRICAL HAZARD: High voltage will be present on the power input module. All service work must be completed by Zimmer Biomet or a qualified biomedical engineer or other person thoroughly familiar with electronic medical devices within the healthcare organization.

TROUBLESHOOTING

SYMPTOM	POSSIBLE CAUSES	SOLUTION
Handpiece will not turn ON and: The battery has blue indicator lights	The handpiece system has detected an error and has stopped all function. Dirty or corroded battery or handpiece terminals. Battery is too hot. The motor has issues. The controller has issues.	Rest the handpiece by disconnecting the handpiece from the power source. Clean the terminals. Allow battery to cool to room temperature. Return for service if the problem persists.
Handpiece will not turn ON and: The battery has no blue indicator lights	The battery is empty. The battery has reached end of life.	Recharge the battery. Replace the battery.
The handpiece is running but the attachment will not move any of the cutting tools	The attachment is not fully inserted into the handpiece. The drivetrain has issues.	Remove and re-insert the attachment. Make sure the attachment is fully inserted and clicks into place. Return for service if the problem persists.
The handpiece is running without any triggers being pressed.	The trigger(s) are contaminated with debris causing sticking. The trigger(s) are collapsed and not returning. The controller has issues.	Clean the triggers per manual cleaning instructions. Add approved lubricant to the trigger shaft(s). Return for service if the problem persists.
The handpiece becomes abnormally hot during use.	The attachment has issues. The motor has issues. The controller has issues.	Replace the suspect attachment. Return for service if the problem persists.
The battery becomes abnormally hot during use.	The battery has issues.	Install battery on battery charger to determine if there are any faults. Replace battery if the problem persists.

SYMPTOM	POSSIBLE CAUSES	SOLUTION
The attachment is stuck and cannot be removed from or installed into handpiece.	The attachment barrel is dirty and contains dried debris. The handpiece attachment socket is dirty and contains dried debris. The attachment has been dropped or damaged. The handpiece has been dropped or damaged.	Clean attachment per the cleaning instructions. Clean the handpiece per the cleaning instructions. Replace the damaged attachments. Return the handpiece for service if the problem persists.
The attachment creates an abnormal noise or vibration when used.	The attachment is misaligned in the handpiece. The attachment has issues.	Remove the attachment and re-install. Replace the damaged attachment.
The handpiece has an abnormal noise or vibration when used.	The handpiece has issues.	Return the handpiece for service if the problem persists.
The handpiece runs but cuts-out under abnormal heavy loading.	The handpiece is experiencing abnormal loading and is entering protection mode. The battery is weak.	Avoid abnormal loading when possible. Install battery on battery charger to determine if there are any faults. Replace battery if the problem persists.
The handpiece slows down under heavy loading.	The handpiece is experiencing extreme loading and reducing the power to avoid entering protection mode. The battery is weak.	Avoid extreme loading when possible. Install battery on battery charger to determine if there are any faults. Replace battery if the problem persists.

EXPECTED SERVICE LIFE INFORMATION

- The *X Series* Power System must be returned every 6 months for inspection and preventive maintenance. Annual factory calibration checks are required to verify continued accuracy. Always follow maintenance requirement.
- The handpiece, attachments and battery charger have an expected service life of 7 years with bi-annual preventive maintenance. Remove product from service at the end of the expected service life. Expected service life based on typical usage of approximately 238 uses a year.
- The battery has an expected service life of 3 years with typical usage. Remove product from service at the end of the expected service life.
- The handpieces have built-in software to protect the internal electronic components. When this protection is activated, the handpiece may not be responsive. To reinitiate the system, disconnect the battery and reconnect it.

RETURN AUTHORIZATION AND REPLACEMENT INFORMATION

- When it is necessary to return the *X Series* Power System for inspection, preventive maintenance or repair within the U.S.A. call 1-800-830-0970 to receive a Return Goods Authorization (RGA) number.
- Outside the U.S.A., contact your local *Zimmer Biomet* representative.
- The instrument must be properly packaged when sent for repair. If the original packaging is no longer available, proper packaging can be requested when the RGA is received.
- A purchase order must accompany all equipment for repair. The customer will be responsible for all shipping charges.

REPLACEMENT PARTS

Please contact [Zimmer Biomet service](#) or your local representative for a list of available replacement parts and order information

To ensure prompt service, please include the following information with your order:

Model Number	Shipping Means (if any)	Shipping Address
Description of Part	Serial Number	
Quantity Desired	Part Number (if known)	

CAUTION: All repairs and preventive maintenance must be done by a *Zimmer Biomet* authorized service center.

WARRANTIES

WARRANTY INSIDE THE U.S.A.

Zimmer Biomet, Inc. warrants the Product (X Series Power System) for one year from date of purchase; excluding batteries which warrants six (6) months. During the warranty period, Zimmer Biomet will repair or replace, at its option, any product which is defective in materials or workmanship or which fails to meet the published specification for that model. This Limited Warranty is made only to the original purchaser of the product and is non-transferable. The remedies described in this Limited Warranty are the exclusive remedies for breach of warranty. **THIS WARRANTY SHALL NOT APPLY TO ANY PRODUCT WHICH HAS BEEN ALTERED OR MODIFIED IN ANY WAY, OR WHICH HAS BEEN SUBJECTED TO MISUSE OR ABUSE.**

DISCLAIMER OF IMPLIED WARRANTIES

The forgoing of Express Limited Warranty is given in lieu of any and all other express or implied warranties. **ZIMMER BIOMET MAKES NO OTHER WARRANTIES INCLUDING THE IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

LIMITATION OF REMEDIES

In no case shall Zimmer Biomet, Inc. be liable for any special, incidental, or consequential damages whether based on breach of warranty or other legal theory whether or not such damages are foreseeable. Some states do not allow limitations on warranties or on remedies for breach in certain transactions. In such states, the limits in this paragraph and the preceding paragraph do not apply.

WARRANTY CLAIMS

In the event of a warranty claim within the warranty period please take the following steps:

- Notify Customer Service Department, Zimmer Biomet Surgical, at 1-800-348-2759 or contact your local Zimmer Biomet representative. Please provide details regarding the nature of the claim and include the product serial number. Upon receipt of this information, Zimmer Biomet will provide a date for service or a return shipping authorization.
- Upon receipt of the shipping authorization, forward the equipment, freight prepaid, to the location specified in the shipping authorization.

Your compliance with these steps will help assure that you receive prompt warranty service for your product.

WARRANTY OUTSIDE U.S.A.

SCOPE OF WARRANTY

Please contact your local Zimmer Biomet Representative for warranty information.

SPECIFICATIONS AND PERFORMANCE

CHARACTERISTICS

~ AC Power Mains Line Voltage Range for Battery Charger	100–240 V ~ (AC), 50/60 Hz. Auto switching.
Input Power – Battery Charger	660 VA
Battery Type	Rechargeable, Li-Ion, 14.4 VDC, 2000 mAh.
Battery Capacity	28.8 Watthours
Battery Recharge Time	1 hour
Power Cord for (Battery Charger)	Type SJT or international equivalent of HO5 VV-F, AWG 16 (1.31 mm ²), 14 ft. (4.27 m).
Power Plug	Hospital grade, 3 prong straight blade, 15 amp or proper equivalent power plug.
Handpiece Battery Line Protection Fuses (Battery Chargers)	Contains no replaceable fuses Contains no replaceable fuses 2 time-delayed 6.3 amp 250 volt fuses
Internal Diagnostics	Program, memory, watchdog timer, overcharge, undercharge, over discharge, motor sensing, battery authentication.
Handpiece Size	Height: 6 in. (12.2 cm) Width: 2.1 in. (5.3 cm) Depth: 5 in. (12.8 cm) Weight: 1.7 lbs (0.77 Kg)
Battery Charger Size	Height: 13.5 in. (34.3 cm) Width: 10.8 in. (27.4 cm) Depth: 8.6 in. (21.8 cm) Weight: 12.4 lbs (5.6 Kg)
Aseptic Battery Housing	Height: 2.4 in. (6.1 cm) Width: 3.4 in. (8.6 cm) Depth: 4.2 in. (10.7 cm) Weight: 0.45 lbs (0.205 Kg)
Aseptic Battery	Height: 1.6 in. (4.1 cm) Width: 3.1 in. (7.9 cm) Depth: 3.5 in. (8.9 cm) Weight: 0.62 lbs (0.280 Kg)
Sterilizable Battery	Height: 2.3 in. (5.9 cm) Width: 3.5 in. (8.8 cm) Depth: 4.1 in. (10.3 cm) Weight: 0.9 lbs (0.390 Kg)

SPECIFICATIONS AND PERFORMANCE

IEC 60601-1 Classification

Class of Protection	Type of protection against electric shock	Handpiece: Internally Powered Equipment Battery charger: Class I
	Degree of protection against electric shock	Type BF applied part.
Handpiece IPX6	Classification according to the degree of protection against ingress of water.	IPX6 – Protection against powerful water jets.
Sterilizable Battery IPX5	Classification according to the degree of protection against ingress of water.	IPX5 – Protection against a sustained, low-pressure water jet spray.
Aseptic Battery IPX1	Classification according to the degree of protection against ingress of water.	IPX1 – Protection against vertically falling water drops
Mode of operation	Intermittent Operation	Forward and Reverse duty cycle of: 1 minute ON / 5 minute OFF; 5 times with a 2 hour rest. Oscillate duty cycle of: 15 seconds ON / 15 seconds OFF, 5 times with a 2 hour rest.
Operator position	Operators position when using handpiece for label legibility.	Maximum of 15.7 in. (40 cm).
Applied part of handpiece	Temperature of the external surface of the handpiece.	Maximum of 118 °F (47.7 °C).

CAUTION: These devices are not suitable for use in the presence of flammable anesthetic or gases.

EMISSIONS/IMMUNITY

The X Series Power System complies with EMC criteria set forth in IEC 60601-1-2. The user of this device should be aware that precautions should be taken in regards to EMC. The device should be installed and used according to the EMC information provided in the instructions for use. See EMC Guidance Tables included in this manual.

WARNING: The reciprocal interference produced by other medical equipment may adversely influence the function of the X Series Power System. When other medical devices are in use, maximize the separation distance and verify that the X Series Power System including the battery charger is working correctly. Do not use in environments with uncontrolled or strong radiofrequency energy, a static magnetic field, or a magnetic field that changes with time. If placed in these environments, the X Series Power System including the battery charger should be observed to verify it is functioning normally.

Cable	Maximum length
~ AC Power Mains power cord	14 ft. (4.27 m)

WARNING: use of an ~ AC Power Mains power cord with a length other than those specified may result in increased emissions and decreased immunity.

ELECTROMAGNETIC COMPATIBILITY (EMC) GUIDANCE TABLES

The following tables provide guidance on needs and installation of the X Series Power System regarding electromagnetic compatibility.

EMC GUIDANCE AND DECLARATION – EM EMISSIONS

Guidance and manufacturer's declaration – electromagnetic emissions.		
The X Series Power System is intended for use in the electromagnetic environments specified below. The customer or the user of the X Series Power System should assure that it is used in such an environment.		
EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
RF emissions CISPR 11 *Only radiated for the handpiece connected to the battery ** Radiated and conducted for the X Series Power System battery charger	Group 1	The X Series Power System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11 *Only radiated for the handpiece connected to the battery ** Radiated and conducted for the X Series Power System battery charger	Class A	The X Series Power System is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2 *Not applicable for the handpiece connected to the battery	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3 *Not applicable for the handpiece connected to the battery	Complies	

EMC GUIDANCE AND DECLARATION – EM IMMUNITY/DISTURBANCES

Guidance and manufacturer's declaration – electromagnetic emissions.			
The X Series Power System is intended for use in the electromagnetic environments specified below. The customer or the user of the X Series Power System should assure that it is used in such an environment.			
IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT – GUIDANCE
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact	± 8 kV contact (230 VAC 50 Hz / battery)	Floors should be wood, concrete or ceramic tile.
	± 15 kV air	± 15 kV air (230 VAC/ battery)	If floors are covered with synthetic materials, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines	± 2 kV for power supply lines (120 VAC 60 Hz / 230 VAC 50 Hz)	~ AC Power Mains power quality should be that of a typical commercial or hospital environment.
*Not applicable for the handpiece connected to the battery			
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	± 1 kV line(s) to line(s) (120 VAC 60 Hz / 230 VAC 50 Hz)	~ AC Power Mains power quality should be that of a typical Commercial or Hospital environment.
*Not applicable for the handpiece connected to the battery	± 2 kV line(s) to earth	± 2 kV line(s) to earth (120 VAC 60 Hz / 230 VAC 50 Hz)	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 *Not applicable for the handpiece connected to the battery	<5% UT (>95% dip in UT) for 0.5 cycle	<5% UT (>95% dip in UT) for 0.5 cycle (120 VAC 60 Hz / 230 VAC 50 Hz)	~ AC Power Mains quality should be that of a typical commercial or hospital environment. If the user of the X Series Power System requires continued operation during ~ AC Power Mains interruptions, it is recommended that the X Series Power System (battery charger or power supply) be powered from an uninterruptible power supply or a battery.
	40% UT (60% dip in UT)	40% UT (60% dip in UT) (120 VAC 60 Hz / 230 VAC 50 Hz)	
	70% UT (30% dip in UT)	70% UT (30% dip in UT) (120 VAC 60 Hz / 230 VAC 50 Hz)	
	<5% UT (>95% dip in UT) for 5 s	<5% UT (>95% dip in UT) for 5 s (120 VAC 60 Hz / 230 VAC 50 Hz)	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U _T is the ~ AC Power Mains voltage prior to application of the test level.			

SPECIFICATIONS AND PERFORMANCE

EMC GUIDANCE AND DECLARATION – EM EMISSIONS/RF

Guidance and manufacturer's declaration – electromagnetic emissions.

The X Series Power System is intended for use in the electromagnetic environments specified below. The customer or the user of the X Series Power System should assure that it is used in such an environment.

Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6 *Not applicable for the handpiece connected to the battery	3 Vrms 150 kHz to 80 MHz	3 V 6 V for ISM Bands	Portable and mobile RF communications equipment should be used no closer to any part of the X Series Power System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance levels in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m	

NOTES:

- At 80 MHz and 800 MHz, the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephone and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the X Series Power System is used exceeds the application RF compliance levels above, the X Series Power System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the X Series Power System.
- b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

EMC GUIDANCE AND DECLARATION – IMMUNITY/SEPARATION DISTANCES

Recommended separation distances between portable and mobile RF communication equipment and the X Series Power System

The X Series Power System is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the X Series Power System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the X Series Power System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in Watts (W)	Separation distance according to frequency transmitter in Meters (m)		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTES:

- At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

ENVIRONMENTAL CONDITIONS

The following are environmental conditions for operation of the X Series Power System:

CONDITION	OPERATION	STORAGE	TRANSPORT
Temperature	50 °F (10 °C)  86 °F (30 °C)	-0.4 °F (-18 °C)  140 °F (60 °C)	-0.4 °F (-18 °C)  140 °F (60 °C)
Relative Humidity	20%  60%	30%  85%	30%  85%
Atmospheric Pressure	800 hPa  1060 hPa	800 hPa  1060 hPa	800 hPa  1060 hPa
Altitude	≤2,000 m	≤2,000 m	--
Pollution Degree	2	2	2

The following are the sounds pressure emission levels and sound power levels according to EU Directive 2006/42/EC:

DEVICE	ATTACHMENT	CUTTING TOOL (Blade)	SOUND PRESSURE LEVEL (LpA) in dBA
Double Trigger Handpiece	Sagittal Saw Attachment	25090127XG1	85.5
Double Trigger Handpiece	Sagittal Saw Attachment	21100147XFW	89.5
Double Trigger Handpiece	Sagittal Saw Attachment	06025060XTF	78.2
Double Trigger Handpiece	Reciprocating Saw Attachment	12090119UR1	90.6
Double Trigger Handpiece	Reciprocating Saw Attachment	12076100UR1	90.3
Double Trigger Handpiece	Reciprocating Saw Attachment	06048119URL	90.2

The following are the vibration emission levels according to EU Directive 2002/44/EC:

DEVICE	ATTACHMENT	CUTTING TOOL (Blade)	HAND TRANSMITTED VIBRATION (m/s ²)	MAXIMUM DAILY EXPOSURE
Double Trigger Handpiece	Sagittal Saw Attachment	25090127XG1	< 2.5	No limitation
Double Trigger Handpiece	Sagittal Saw Attachment	21100147XFW	< 2.5	No limitation
Double Trigger Handpiece	Sagittal Saw Attachment	06025060XTF	< 2.5	No limitation
Double Trigger Handpiece	Reciprocating Saw Attachment	12090119UR1	3.97	3 hours 10 minutes
Double Trigger Handpiece	Reciprocating Saw Attachment	12076100UR1	5.7	1 hour 32 minutes
Double Trigger Handpiece	Reciprocating Saw Attachment	06048119URL	< 2.5	No limitation

NOTES:

- The double trigger device is used for the sound pressure and vibration testing. There is no difference between the double and single trigger devices. The worst case cutting tool was utilized for the test to report maximum levels.
- The maximum daily exposure to hand vibration is considered for a cumulative time during a 24 hours period.

SYMBOLLOGY

Below are graphical Warning, Caution and other Symbols indicated on the Power System and within this manual.

Graphic	What it means
	Signifies to follow instructions for use (IFU).
	Consult instructions for use.
	Do not immerse.
	Fuse.
	Medical - General Medical Equipment as to Electrical Shock, Fire and Mechanical hazard only, in Accordance with ANSI/AAMI ES60601-1:2005+AMD1:2012 and CAN/CSA-C22.2 No.60601-1:2008
	This product contains one or more toxic or hazardous substances or elements. The Environmental Protection Use Period on the logo refers to the period in years which toxic or hazardous substances or elements contained in the product will not leak or mutate under normal operating conditions.
	This product contains electrical or electronic materials. The presence of these materials may, if not disposed of properly, have potential adverse effects on the environment. Presence of this label on the product means it must not be disposed of in normal household waste and must be disposed of separately. To find out how to properly dispose of this product, please contact your local Zimmer Biomet Representative.
	This product contains a Li-ion battery that must be recycled.
	CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.
	Conform to the European Medical Device Directive 93/42/CEE + 2007/47/CE for class IIa medical devices (notified body number).
	Conform to the European Medical Device Directive 93/42/CEE + 2007/47/CE for class I medical devices.
	Type BF applied part.
	Identifies any terminal which is intended for connection to an external conductor for protection against electric shock in case of a fault, or the terminal of a protective earth (ground) electrode.
	Potential equalization connection. Identifies the terminals which, when connected together, bring the various parts of an equipment or of a system to the same potential, not necessarily being the earth (ground) potential, e.g. for local bonding.
	Atmospheric pressure limitation.
	Temperature limit.
	Humidity limitation.
	Indicates hazards arising from dangerous voltages.
	Signifies a general warning.
	Signifies a general caution.

Graphic	What it means
	Used to identify the devices sold in nonsterile conditions.
	Indicates authorized representative (REP) in the European Community (EC).
	Indicates importer in the European Community (EC).
	Indicates catalog, reorder or reference number.
	Indicates equipment serial number.
	Indicates equipment lot number.
	Indicates Unique Device Identifier number.
	Quantity contained in package.
	Indicates a product intended to be used as a medical device.
 Note: Can be used in combination with date of manufacture.	Indicates manufacturer and is accompanied by the name and address of the manufacturer.
 YYYY-MM	Indicates date of manufacture and is accompanied by a date.
	Indicates suitability for alternating current only.
	Indicates suitability for direct current only.
	Indicates equipment or systems in the medical electrical area that include RF transmitters or that intentionally apply RF electromagnetic energy for diagnosis or treatment.

TEST AND EVALUATION REPORT

Testing Organization: _____		Test before putting into service (reference value) ☐	
Name of Testing Person: _____		Recurrent test ☐	
Responsible Organization: _____		Test after repair ☐	
Equipment: _____		ID – Number: _____	
Type: _____		Serial Number: _____	
Manufacturer: Zimmer Biomet		Class of Protection: Class I ☐ Class II ☐ Battery ☐	
Applied Part Type: 0 ☐ B ☐ BF ☐ CF ☐	Mains Connection: Permanently Installed ☐ Non-detachable Power Cord ☐ Detachable Power Cord ☐		
Accessories: _____			
Measuring Equipment: _____			
Test			Results
Visual Inspection			Pass ☐ Fail ☐
Measurements:	Specification	Measured	
Protective earth resistance (cord attached):	≤300 mΩ	_____ mΩ	Pass ☐ Fail ☐
Equipment Leakage:	< 10mA	_____ mA	Pass ☐ Fail ☐
Insulation Resistance – Battery Charger:	>70MΩ	_____ MΩ	Pass ☐ Fail ☐ N/A ☐
Insulation Resistance – Handpiece:	>70MΩ	_____ MΩ	Pass ☐ Fail ☐ N/A ☐
Functional Test:	Specification		
Mode switch at OFF – Upper trigger	Trigger locked and does not turn ON		Pass ☐ Fail ☐
Mode switch at OFF – Lower trigger	Trigger locked and does not turn ON		Pass ☐ Fail ☐ N/A ☐
Mode switch LEFT – Upper trigger	Trigger unlocked and motor turns ON		Pass ☐ Fail ☐
Mode switch LEFT – Lower trigger	Trigger unlocked and motor turns ON		Pass ☐ Fail ☐
Releasing trigger(s) causes motor to stop	Motor stops when trigger(s) released		Pass ☐ Fail ☐
Mode switch RIGHT – Upper trigger	Trigger unlocked and motor turns ON		Pass ☐ Fail ☐
Mode switch RIGHT – Lower trigger	Trigger unlocked and motor turns ON		Pass ☐ Fail ☐ N/A ☐
Releasing trigger(s) causes motor to stop	Motor stops when trigger(s) released		Pass ☐ Fail ☐
Deficiency / Notes: _____			
Overall Assessment: No safety or functional deficiencies were detected ☐ No direct risk, deficiencies detected may be corrected on short term ☐ Equipment shall be taken out of service until deficiencies are corrected ☐			
Next recurrent test necessary in:		1 month ☐ 3 months ☐ 6 months ☐ 12 months ☐	
Name: _____ Signature: _____ Date: _____			

CONTACT INFORMATION

To obtain part or additional information regarding your unit, contact your Zimmer Biomet distributor or write or phone:

Mail	Telephone
Zimmer Surgical SA 3 Chemin du Pré-Fleuri 1228 Plan-les-Ouates Switzerland	• Inside the U.S.A., call the Zimmer Biomet Biomet Customer Service Department at 1-800-348-2759. • Outside the U.S.A., contact your local Zimmer Biomet sales representative.

REPORTING PROBLEMS

The user should report any suspected serious incident related to the X Series Power System by informing the manufacturer and the competent authority of the state or the country in which the serious incident has occurred.

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©2020 Zimmer Biomet